


# SLIMPLICITY® ANTERIOR CERVICAL PLATING SYSTEM



**Low Profile Plate (2mm)**  
designed to minimize tissue disruption and post-op discomfort

**Pre-Contoured**  
to address patient anatomy



**Central Graft Windows**  
facilitate unimpeded graft visualization

**One Step Locking**  
facilitates secure application

## IMPLANT SIZING

### Plate Sizing

Lengths\*

- 1 Level - 7 sizes, 20-32mm (2mm increments)
- 2 Level - 8 sizes, 37-55mm (3mm increments)
- 3 Level - 8 sizes, 54-75mm (3mm increments)
- 4 level - 6 sizes, 69-89mm (4mm increments)


\*Measured end-to-end

Width 17mm

Waist 14mm


### Fixed & Variable Screw Sizes

- 4.0 and 4.5mm Diameters
- 12-18mm Lengths  
(2mm increments)



Single Tray Layout minimizes case complexity

## PRECISION SPINE SHURFIT® CERVICAL INTERBODY CAGES



- **Large Contact Area**  
designed to optimize vertebral body support and minimize risk of subsidence
- **Large Graft Area**  
enables excellent biological coverage
- **Trapezoidal Design**  
allows for proper anterior placement
- **Tantalum Markers**  
strategically placed to aid interbody positioning

### SIZING

- 14mm Wide x 12mm Deep
- 16mm Wide x 14mm Deep
- Heights - 5-12mm (1mm increments)
- 5° & 10° Lordotic Options

**Precision Spine, Inc.**  
2050 Executive Drive, Pearl, MS 39208  
Customer Service: 1.888.241.4773  
Phone: 601.420.4244  
Toll Free: 877.780.4370  
Fax: 601.420.5501  
www.precisionspineinc.com

**PRECISION SPINE®**  
CE 2797